



ES-150V Pro Electrosurgical system

Designed to meet the highest veterinary standards, the ES-150V Pro delivers precise, safe, and dependable performance.

Main Feature

- Full-Glass Design
- Color touchscreen operation panel — flexible and easy to use
- 7 Working modes (5 monopolar and 2 bipolar)
- Large blood vessel sealing function
- Electrosurgical pens and foot switch control
- Recent mode, power, and other parameter display
- Quick power adjustment
- Intermittent cut and coagulation modes





Mode			Max Output Power(W)	Load Impedance(Ω)	Modulation Frequency (kHz)	Max Output Voltage (V)	Crest Factor
Monopolar	Cut	Pure Cut	120	500	--	1300	1.8
		Blend 1	100	500	20	1400	2.0
		Blend 2	100	500	20	1300	2.0
	Coag	Spray	90	500	12-24	4800	6.3
		Forced	60	500	25	4800	6.2
Bipolar	Vessel Sealing		--	100	20	700	1.9
	Standard		60	100	20	700	1.9

Ordering Information

Model	Description
	Electrosurgical Unit
	Monopolar Foot-Switch
	Bipolar Foot-Switch
	Hand-Switch Pencil, Disposable
	Patient Return Electrode Rods (10mm) With Cable, Reusable
	Maryland Jaw Open and Laparoscopic Sealer, 5mm-37cm



T. +1 949-679-2881
F. +1 949-679-2882



www.medicattechvet.com



Support@medicatechvet.com