



ES-200V Electrosurgical Generator

Multifunctional electrosurgical generator with a wide range of application departments and a very high cost performance

Main Feature

- 3 Monopolar Cut Modes: Pure Cut, Blend 1, Blend 2
- 3 Coagulation Modes: Spray Coag, Forced Coag, Soft Coag
- 2 Bipolar Outputs: Standard, Precise
- CQM (Contact Quality Monitoring) System
- Electrosurgical pens and foot switch control
- Recent mode, power, and other parameter display
- Quick adjustment of power and volume
- Intermittent cut and coagulation modes





Mode			Max Output Power(W)	Load Impedance(Ω)	Modulation Frequency (kHz)	Max Output Voltage (V)	Crest Factor
Monopolar	Cut	Pure Cut	200	500	--	1300	1.8
		Blend 1	200	500	20	1400	2.0
		Blend 2	150	500	20	1300	2.0
	Coag	Spray	120	500	12-24	4800	6.3
		Forced	120	500	25	4800	6.2
		Soft	120	500	20	1000	2.0
Bipolar	Standard		100	100	20	700	1.9
	Precise		50	100	20	400	1.9

Ordering Information

Model	Description
	Electrosurgical Unit
	Monopolar Foot-Switch
	Bipolar Foot-Switch
	Hand-Switch Pencil, Disposable
	Plastic and Aesthetic Surgery / Dermatology / Oral / Maxillofacial Surgery
	Patient Return Electrode Without Cable, Split, for Adult, Disposable
	Connecting Cable for Patient Return Electrode(Split) , 3m, Reusable



T. +1 949-679-2881
F. +1 949-679-2882



www.medicattechvet.com



Support@medicatechvet.com