



Medicatech VET[®]
A Medicatech Company

Apogee 1000 Lite Veterinary Digital Color Doppler Ultrasound Imaging System





Medicatech VET[®]

A Medicatech Company

PORTABILITY, FUNCTIONALITY, AFFORDABILITY

Apogee 1000V is a laptop-designed veterinary ultrasound system that provides you with an all-in-one veterinary solution. Incorporating outstanding image quality, multifunctional veterinary packages, and advanced transducers, the Apogee 1000V Lite enhances general practice capabilities while supporting your financial goals.



COMPREHENSIVE VETERINARY SOLUTIONS

Apogee 1000V Lite offers complete veterinary measurement packages with multiple presets covering abdomen, cardiology, reproduction, small parts, tendonS, etc.

SCANNING MODES : B, Color, PW, CW, M, AMM, 3D Mode, etc.

VETERINARY EXAM PACKAGES : Measurement, report and body mark dedicated for canine, feline, equine, bovine, ovine and more.

AUTOMATIC MEASUREMENT : Auto Trace, Auto Volume, Semi Auto EF, etc.

COMPETENT EXAM TOOLS : TDI, Color M, Elastography, Needle Enhancement, Panoscope, VS Flow, etc.





Medicatech VET[®]

A Medicatech Company

HIGH PERFORMANCE TRANSDUCERS

A wide range of transducers is available to support comprehensive whole-body examinations for companion animals, equine, and farm animals.

For Companion Animals



L8
Linear Probe



C6
Micro Convex Probe



L1FC
Backfat Probe



L8LVN
Veterinary Transrectal Probe

For Equine



C3
Convex Probe



P3
Phased Array Probe




L5LV
Veterinary Trans-rectal Probe




KIDNEY



SPLEEN



INTESTINE



HEART



T. +1 949-679-2881
F. +1 949-679-2882



www.medicattechvet.com



Support@medicatechvet.com