



MEDICATECH VET®

A Medicatech Company

Veterinary Reusable Flexible Endoscope

Veterinary Single-use Flexible Endoscope

Veterinary Single-use Gastroscope





Ordering information

Veterinary Reusable Flexible EndoScope

Product	Item code
Veterinary reusable flexible endoscope	
VT-Q28	2402B-PA001
VT-Q30	2402B-PA002
VT-Q40	2402B-PA003
VT-Q52	2402B-PA004
VT-Q58	2402B-PA005
VT-Q28S	2402B-PA006

Veterinary Single-use Flexible EndoScope

Product	Item code
Veterinary single-use flexible endoscope	
VBR-M22	2016B-PA055
VBR-M25	2016B-PA061
VBR-M32	2016B-PA054
VBR-M40	2016B-PA053
VBR-M50	2016B-PA052
VBR-M58	2016B-PA051
VBR-M40H	2016B-PA050
VBR-M52H	2016B-PA049
VBR-M58H	2016B-PA048
VBR-M50-10	2016B-PA060
VBR-M58-15	2016B-PA058
VBR-M52H-10	2016B-PA059
VBR-M58H-15	2016B-PA057


Medical Image Processor

Product	Item code
Medical Image processor	
MS-8	2288B-PA001
Medical Image processor	
VL35-V	200-0410-00




SCOPE OF APPLICATION SMALL ANIMALS

VT-Q28/VT-Q28S/VT-Q30 ($\Phi 2.8\text{mm}/\Phi 3.0\text{mm}$)




Feline

VT-Q40/VT-Q52 ($\Phi 4.0\text{mm}/\Phi 5.2\text{mm}$)



Canine

VT-Q28/VT-Q28S/VT-Q30/VT-Q40/VT-Q52, OD(2.8mm-5.2mm), ID(1.2mm-2.6mm) is suitable for small and medium animals' bronchoscopy, rhinoscopy, oesophagus, foreign body retrieval, gastroscopy, cystoscopy, vaginoscopy, etc.



Equine

with its long working length and big working channel, VT-Q58 is ideal for equine's upper and lower airway endoscopy, sinuscopy, vaginoscopy and hysteroscopy, etc



MEDICATECH VET®
A Medicatech Company



One-stop management

1

**Veterinary
Reusable Flexible Endoscope**

2

**Veterinary
Single-use Flexible Endoscope**

3

**Veterinary
Single-use Gastroscope**



MEDICATECH VET®
A Medicatech Company

Reusable Flexible Endoscope



Display

- 15.6 inch touch screen
- All in one system



120° Left/Right Rotation

- Relieve fatigue of wrist rotation



Distal tip part

- Standard 160K pixel



Function

- Multi-function Buttons: Photo & video recording
- Suction Button
- Working channel



Model	O.D.	I.D.	Bending Angle(up/down) Rotate (left/right)	Field of View	Working Length
VT-Q28	2.8mm	/	210° Up / 130° down, 120° Left / right	110°±10%	650mm
VT-Q28S	2.8mm	1.2mm	210° Up / 130° down, 120° Left / right	110°±10%	650mm
VT-Q30	3.0mm	1.2mm	180° Up / 130° down, 120° Left / right	110°±10%	650mm
VT-Q40	4.0mm	2.0mm	180° Up / 130° down, 120° Left / right	110°±10%	1000mm
VT-Q52	5.2mm	2.6mm	180° Up / 130° down, 120° Left / right	110°±10%	1000mm
VT-Q58	5.8mm	3.0mm	180° Up / 130° down, 120° Left / right	110°±10%	1500mm



MEDICATECH VET®
A Medicatech Company

Single-use Flexible Endoscope INTENT/ MEGA

Multiple model selection



Chip on tip
CMOS+LED



Intensive Bending
210° up/down



Lightweight
≤0.28 Kg



Rotation
120° left/right

Type	Model	O.D.	I.D.	Bending Angle(up/down) Rotate (left/right)	Field of View	Working Length
INTENT	VBR-M22	2.2mm	/	210° Up/down, 120° Left/right	120°±10%	600mm
	VBR-M25	2.5mm	1.2mm	210° Up/down, 120° Left/right	120°±10%	600mm
	VBR-M32	3.1mm	1.2mm	210° Up/down, 120° Left/right	120°±10%	600mm
	VBR-M40	4.0mm	2.0mm	210° Up/down, 120° Left/right	120°±10%	600mm
	VBR-M50	5.0mm	2.8mm	210° Up/down, 120° Left/right	120°±10%	600mm
	VBR-M58	5.8mm	3.4mm	210° Up/down, 120° Left/right	120°±10%	600mm
	VBR-M50-10	5.0mm	2.8mm	210° Up/down, 120° Left/right	120°±10%	1000mm
	VBR-M58-15	5.8mm	3.4mm	210° Up/down, 120° Left/right	120°±10%	1500mm
MEGA	VBR-M40H	4.0mm	1.2mm	210° Up/down, 120° Left/right	120°±10%	600mm
	VBR-M52H	5.0mm	2.2mm	210° Up/down, 120° Left/right	120°±10%	600mm
	VBR-M58H	6.0mm	2.8mm	210° Up/down, 120° Left/right	120°±10%	600mm
	VBR-M52H-10	5.0mm	2.2mm	210° Up/down, 120° Left/right	120°±10%	1000mm
	VBR-M58H-15	6.0mm	2.8mm	210° Up/down, 120° Left/right	120°±10%	1500mm

View clearer, do better.

MS-8 15.6 inches touch screen Image Processor

The all-around upgrade digital video processor for HugeMed adopts a new operating system, which functions as a central system to facilitate easier reading and operating in various endoscopy applications. Especially important for occasions where mobility is needed in hospitals.


Touch
Control menu

Patient information
Management



All-in-One system

Scope Operating
Management




15.6 Inch



1920*1080 Full HD



Touch Screen



170° View Angle

DICOM & WiFi

Utilizing DICOM protocol, information can seamlessly transfer to the hospital PACS system, effectively enhancing inter-departmental collaboration efficiency.

Image adjustment function

LED level/sharpness/color tone/contrast provide you with better visual experience.

Patient information management


Operator can edit medical records then generate a graphic report.

High compatibility

Can support the connections of HugeMed single-use flexible endoscope and reusable flexible endoscope.

Third-party LCD monitor connection

HDMI/DVI/SDI image output ability.





MEDICATECH VET®
A Medicatech Company

Gastroscope System

Single-use gastroscope



■ Specification


O.D.	I.D.	Bending Angle(up/down) Rotate (left/right)	Working Length
6.0mm	2.8mm	210° Up/down, 140° Left/right	1200mm
7.5mm	2.8mm	210° Up/down, 140° Left/right	1500mm
10.0mm	3.8mm	210° Up/down, 140° Left/right	1500mm
10.0mm	3.8mm	210° Up/down, 140° Left/right	3200mm



Single-use
Sterilized



Ergonomic Design
Ease-of-use



Chip on tip
CMOS+LED



T. +1 949-679-2881
F. +1 949-679-2882



www.medicattechvet.com



Support@medicatechvet.com