



MEDICATECH VET®
A Medicatech Company

Hand-Held Dental Generators

MasterDent Pro





► System Specifications

X-ray generator

Tube voltage

Tube current

Tube type

Tube focal spot

Target angle

Exposure time (increase 0.05sec)

Power requirement

Weight

Power supply

Inverter type

50 ~80 kV (adjustable)

1 ~ 5 mA (adjustable)

Stationary

0.8mm

16°

0.05~1.35sec


DC 15.0 ~ 16.8V

3.6kg

Rechargeable battery (16.8V) -Inner
battery pack –External battery pack (At
the same time charging and using)

Display

4.8" TFT-LCD(Radiation image display)





T. +1 949-679-2881
F. +1 949-679-2882



www.medicattechvet.com



Support@medicatechvet.com