



**Medicatech VET<sup>®</sup>**  
A Medicatech Company

# Veterinary Mobile Dental System





## Technical Parameters

### Product size

44\*35\*72cm

### Machine parameters and description

Power: 0.6KW  
Voltage: 220V/50HZ  
Exhaust volume:  $\geq 118$ L/min  
Noise value:  $48 \pm 2$  dB(A)  
Operating pressure: 0-6 Bar  
Water bottle/sewage bottle Capacity: 650mL/1000mL  
Mobile phone tube :2 (2-hole / 4-hole optional)  
Silent air compressor :1  
Three guns: one  
Saliva inhalation :1  
curing Light :1  
Tooth cleaning machine :1  
Stainless steel gas storage cylinder :1  
Clean water bottle :1  
Collection bottle :1  
Product size: 44\*35\*72cm  
Packing size: 53\*45\*92cm(including tray)  
Net/gross weight: 35/45kg  
Optional: High speed or low speed phone





T. +1 949-679-2881  
F. +1 949-679-2882



[www.medicattechvet.com](http://www.medicattechvet.com)



[Support@medicatechvet.com](mailto:Support@medicatechvet.com)