



Medicatech VET®
A Medicatech Company



Stainless Steel Wet disposal Treatment Table





Technical Parameters

Material	whole 304 Stainless
Features	<ul style="list-style-type: none">• L1500XW600XH880mm(59"Lx23"Wx34.5" H)• The surface of the disposal platform is made of all round corners• without any blind angle.• The cleaning is convenient with the stainless steel top. The water• drain is made of stainless steel filter wool to assure water will not be• blocked.• Equipped with stainless steel retractable water tap, 2 water outlet• modes.• Thick adjustable stainless steel sliding guide rail drawer• Includes one large storage cabinet, and three drawers.• Base & Table surface: 304 stainless steel• Grid is made of 8mm• diameter main frame, and 4mm• diameter horizontal rod, cross spot• welding.





T. +1 949-679-2881
F. +1 949-679-2882



www.medicattechvet.com



Support@medicatechvet.com