



Medicatech VET[®]
A Medicatech Company

Veterinary Motorized Elevating/ Tilting Surgical Table





Medicatech VET[®]

A Medicatech Company



Technical Parameters

Material

whole 304 Stainless

Features

- Electric lifting
- Standard with Detachable instrument tray and IV pole
- The height of the table is controlled by electric, 760–1060 mm, and there is a movable sink in the middle
- The right-and-left tilt is -5° / 15° , by knob manual control
- The fore-and-aft tilt is 45° , by automatic control
- Both side have the lock shelf
- The table board is combined by two boards, the middle of table board has a movable sink
- Mobile castor
- Tabletop temperature can be adjusted $0-60^{\circ}\text{C}$, automatic constant temperature function, can adjust the required temperature at will. Built-in high temperature protector, safe and reliable
- Equipped with binding belt, simple and firm.





T. +1 949-679-2881
F. +1 949-679-2882



www.medicattechvet.com



Support@medicatechvet.com