



Medicatech VET[®]
A Medicatech Company



Wood Frame Wet disposal Treatment Table






Technical Parameters

Features

- L1500XW600XH880mm(59"Lx23"Wx34.5" H)
- The surface of the disposal platform is made of all round corners
- without any blind angle.
- The cleaning is convenient with the stainless steel top. The water
- drain is made of stainless steel filter wool to assure water will not be
- blocked.
- Equipped with stainless steel retractable water tap, 2 water outlet
- modes.
- Thick adjustable Stainless steel sliding guide rail drawer
- Includes one large storage cabinet, and three drawers.
- Base : Wood, Table surface: 304 stainless steel Grid is made of
- 8mm diameter main frame, and 4mm diameter horizontal rod, cross
- spot welding.





T. +1 949-679-2881
F. +1 949-679-2882



www.medicattechvet.com



Support@medicatechvet.com